

VirusClear VF

Virus Removal Filter





◀ Viruclear PVF Virus Removal Capsule Filter With PVDF Shell

Viruclear VF

Reliable & Effective Virus Removal Filter

At present, in the common drug virus removal process, virus removal filtration can be used as an effective physical step for virus removal, due to its advantages such as simple and gentle operation, clear mechanism, easy validation, it has been widely used for virus safety assurance in purification processes.

At present, the most widely used is the virus removal filter with a cut-off pore size of 20nm.

Virus retention is achieved by size exclusion mechanism, and its mechanism of action is based on virus larger than the pore size of the filter membrane, the product passes through the filter membrane and the virus is retained, which can be removed viruses with an envelope or without envelopes, and at the same time based on the molecular size retention mechanism, in addition to virus filters, the 20nm interception pore size can not only be used to effectively remove parvoviruses, but also it can also effectively remove mouse leukemia virus with a diameter of 80-100nm and other types of virus particles

Cobetter Viruclear VF adopts a continuous gradient structure of PES membrane, which greatly increases the foulant holding space, and has a strong retention capacity for large-particle impurities, e.g., aggregates. The downstream side is composed of a 20nm-resolution filter membrane with a uniform pore size distribution, which can achieve a robust retention capacity for parvovirus and better response to extreme filtration test challenges like pressure interruption. Viruclear VF virus removal filter membrane adopts a highly hydrophilic formula, which effectively improves the hydrophilicity of the PES membrane, and has a higher flux and capacity on the filtration for high-concentration, more hydrophobic materials. At the same time, Viruclear VF virus removal filter can be used for robust linear process scale-up, and online integrity testing. Cobetter has independent intellectual property rights and has applied for many inventions.

Viruclear VF virus removal filter membrane adopts a high hydrophilic formula, which effectively improves the hydrophilicity of PES membrane, and has a higher flux and capacity. At the same time, Viruclear® VF virus removal filters enable robust linear process scale-up and can be tested for integrity in-line.



▲ Viruclear VF Virus Removal Membrane Cassettes

Features

- PES membrane has a unique surface hydrophilic improvement, which endows virus removal membranes with high filtration capacity and high recovery
- Asymmetric membrane structure, high porosity, which endows the membrane with high flux
- Water flux greater than 600LMH@30psi 25°C
- Narrow pore size distribution gives the membrane good separation performance
- Robust parvovirus removal capacity
- Mode virus removal rate $\geq 4.0\log$
- Maintain robust virus removal capacity even under process interruption
- Good chemical stability
- Easy to install, use and test

Typical Application

Monoclonal antibody

Bispecific monoclonal antibody


Fc fusion protein

Nanobody


Recombinant Covid-19 vaccine

Small molecule recombinant proteins (< 180 kD), etc.

Virus Removal Filter Membrane Pore Size Distribution



Viruclear VF membrane pore size distribution



SEM image of virus removal filter membrane

Product Specifications

	Product Specifications	Effective Filtration Area	Application
Application	DS series	2.5/2.8cm ²	For process development and virus clearance validation research
	Membrane filter	4.1cm ²	For process development and virus clearance validation research
	Pilot device	0.017/0.07/0.22m ²	≤200L
	Production device	0.50/1.50m ²	>200L
Nylon Pre-filtration	DS series	3.4cm ²	For process development and virus clearance validation research
	Pilot device	0.025/0.12/0.30m ²	≤200L
	Production device	0.60/1.20/1.8m ²	>200L
Depth Pre-filtration	DS series	4.5cm ²	For process development and virus clearance validation research
	Pilot device	0.027/0.15/0.4m ²	≤200L
	Production device	0.92/1.1m ²	>200L


Material of Products

	Specifications	Material of Membrane	Format	Shell	Effective Filtration Area	Accessories Material
Removal Virus	DS series	PES	Syringe filter	Polypropylene	2.5/2.8cm ²	-
	Pilot device	PES	Disposable capsule filter	Polyvinylidene fluoride	0.017/0.07/0.22m ²	Silicone, Polypropylene
			Membrane Cassette	Silicone	0.018/0.08/0.26m ²	Polypropylene
	Production device	PES	Disposable capsule filter	Polyvinylidene fluoride	0.50/1.50m ²	Silicone, Polypropylene
Membrane Cassette			Silicone	0.50/1.50m ²	Polypropylene	
Nylon Pre-filtration	DS series	Nylon	Syringe filter	Polypropylene	3.4cm ²	Silicone
	Pilot device	Nylon	Disposable capsule filter	Polypropylene	0.025/0.12/0.30m ²	Silicone
	Production device	Nylon	Disposable capsule filter	Polypropylene	0.60/1.20/1.80m ²	Silicone
Depth Pre-filtration	DS series	Cellulose, diatomite	Syringe filter	Polypropylene	4.5cm ²	Silicone
	Pilot device	Cellulose, diatomite	Disposable capsule filter	Polypropylene	0.027/0.15/0.4m ²	Silicone
	Production device	Cellulose, diatomite	Disposable capsule filter	Polypropylene	0.92/1.1m ²	Silicone


Virus Clearance Filtration Process Development Service

Our technical engineers will work with customers to optimize the process parameters during the virus clearance filtration process to obtain a robust, efficient and economical filtration operation.


Virus Removal Filter Filterability Testing Device



Virus Removal Filter Assembly Sketch Map




Optimization Strategy of Virus Clearance Filtration Process



Customer Application Case

Test conditions: mAb concentration 12g/L; constant pressure 30psi; series filtration with PDF depth filter; After filtration, chasing with 30L/m² Buffer



Filtration performance	Viruclear VF
Filtration time (min)	120
Flux decay at filtration endpoint(%)	27
Capacity at filtration endpoint L/m ²	775
Average flux(LMH)	380
Recover(%)	99

Order Information

Viruclear PDT Virus Removal Pre-Filter (Depth Filter)

V P

Application

VP Viruclear PDF



D T

Material of Membrane

DT Depth filter

Type

DS Syringe filter
(4.5cm²)

Quantity/Package

N9 Only for syringe filter,
9pcs / pk
N1 1pcs / pk

P

Market

P Biopharmaceutical

Viruclear PDF Virus Removal Pre-Filter (Depth Filter)

V P

Application

VP Viruclear PDF



D F

Material of Membrane

DF Depth filter

Type

L4 L08TT(0.027m²)
SA CSCD(0.15m²)
SC CSCE(0.4m²)

Quantity/Package

N1 1pcs / pk

P

Market

P Biopharmaceutical

Viruclear PDF Virus Removal Pre-Filter (Depth Filter)

V P

Application

VP Viruclear PDF



D F

Material of Membrane

DF Depth filter

Type

SB CSCB(0.92m²)
SM CSCM(1.1m²)

Quantity/Package

N1 1pcs / pk

P

Market

P Biopharmaceutical

Order Information

Viruclear PNY Virus Removal Pre-Filter (Nylon Filter)

V P

Application

VP Viruclear PNY

N Y

Material of Membrane

NY mNylon,
Modified nylon
membrane

Type

DS Syringe filter
(3.4cm²)

Quantity/Package

N9 Only for syringe filter,
9pcs / pk
N1 1pcs / pk

P

Market

P Biopharmaceutical



Viruclear PNY Virus Removal Pre-Filter (Nylon Filter)

V P

Application

VP Viruclear PNY

N Y

Material of Membrane

NY mNylon,
Modified nylon
membrane

Type

C02 250cm²
L02 0.12m²
L05 0.30m²

Type

TT 3/4" TC

Quantity/Package

N1 1pcs / pk

P

Market

P Biopharmaceutical



Viruclear PNY Virus Removal Pre-Filter (Nylon Filter)

V P

Application

VP Viruclear PNY

N Y

Material of Membrane

NY mNylon,
Modified nylon
membrane

Type

L10 0.60m²
L20 1.20m²
L30 1.80m²

Type

SS 1 1/2" TC

Quantity/Package

N1 1pcs / pk

P

Market

P Biopharmaceutical



Order Information

Viruclear VF Virus Removal syringe filter with PP shell

V F

Application

VF Virus filtration



M P

Material of Membrane

MP mPES,
Hydrophilic modified
PES membrane

Type

DS Syringe filter (2.5cm²)

DB Syringe filter (2.8cm²)

Quantity/Package

N9 Only for syringe filter,
9pcs/pk (3 batches,
3 pcs per batch)

N1 1pcs / pk

P

Market

P Biopharmaceutical

Viruclear VF Virus Removal membrane filter

V F

Application

VF Virus filtration



M P

Material of Membrane

MP mPES,
Hydrophilic modified
PES membrane

Type

FS Circular membrane,
should be installed in
holder to use

Membrane Diameter

25 Only for membrane filter,
diameter 25mm

P

Market

P Biopharmaceutical

Viruclear PVF Virus Removal capsule filter with PVDF shell

V F

Application

VF Virus filtration

E

Material of
Membrane

E PES

S

Configuration
Code

S Standard

Type

LA Lab(0.017m²&0.07m²)

FL Flow(0.22m²0.50m²& 1.50m²)

Filtration Area

002 0.017m²

008 0.07m²

026 0.22m²

050 0.50m²

150 1.50m²

Quantity/Package

N1 1pcs / pk

P

Market

P Biopharmaceutical



Viruclear VF Virus Removal membrane cassettes

V F

Application

VF Virus filtration

M P

Material of Membrane

MP mPES,
Hydrophilic modified
PES membrane

Type

A1 Pilot device 1
(0.018m²)

A2 Pilot device 2
(0.08m²)

B1 Amplification device 1
(0.26m²)

B2 Amplification device 2
(0.5m²)

B3 Amplification device 3
(1.5m²)

Quantity/Package

N1 1pcs / pk

P


Market

P Biopharmaceutical



Our Mission

Through Excellent Products & Sustainable Innovative Solutions,
We Help Customers Solve Process Problems & Increase Yield.



Please contact us for more information

Hangzhou Cobetter Filtration Equipment Co.,Ltd.

Sales Add 16 F, Building 1, Zicheng International Innovation Center,
No.39 Jincheng Road, Xiaoshan District, Hangzhou 311215, China

Factory Cobetter Park, Heshang New Material Industrial Park,
Xiaoshan District, Hangzhou 311265, China

Tel +86 400-070-4266

Fax +86-571-87704256

www.cobetterfiltration.com